



Sterile Sharps Bins

A range of sharps bins, individually double wrapped and sterilised for your convenience. The bins are clean, sterile and ready for use eliminating the need for excessive spraying and time consuming wiping. All of the bins have needle and vacutainer removers for the safe removal and disposal and transportation of sharps.

0.6 Litre Sharps Bin

- ✓ Conveniently fits lead containers in radio pharmacy applications

2.5 & 5 Litre Sharps Bins

- ✓ Yellow lid for medicinal waste
- ✓ Standard and large sizes

Cyto Sharps Bins

- ✓ Purple lid and labelling to comply with EU directives for the safe handling of cytotoxic waste

Order code	Description	Size	Sterile	Sales Unit
XX0478	Sharps Bin	0.6 Litre	Sterile	12 pcs/pk
XX0472	Sharps Bin	2.5 Litre	Sterile	12 pcs/pk
XX0471	Sharps Bin	5 Litre	Sterile	8 pcs/pk
XX0605	Cyto Sharps Bin	5 Litre	Sterile	8 pcs/pk
XX0478C	Cyto Sharps Bin	0.6 Litre	Sterile	12 pcs/pk
XX0472C	Cyto Sharps Bin	2.5 Litre	Sterile	12 pcs/pk

For more information and to place an order contact our expert sales team free on

0800 0328 428

or visit our website www.helapet.co.uk