



Berner Gloves

Berner double dipped latex gloves offer excellent protection against cytotoxic preparations and are extra strong to resist tearing on the collars of isolators and cabinets.

Extensive test and breakthrough data is available for the cytotoxics most commonly used in aseptic units and oncology. The gloves are hand specific with a beaded cuff for comfort and security, and are supplied sterile and aseptically presented in pairs in a cleanroom pouch.

- ✓ **Double dipped latex to give longer breakthrough times for cytotoxics**
- ✓ **Hand specific for comfort**
- ✓ **Beaded cuff to resist tearing on collars of isolators**
- ✓ **Sterile and aseptically presented**
- ✓ **Clean peel pouch for pharmaceutical Cleanrooms**
- ✓ **CE marked**

Order Code	Description	Size	Sterile	Sales Unit
GL3010	Berner Glove	6.5	sterile	50pcs/pack
GL3015	Berner Glove	7.0	sterile	50pcs/pack
GL3020	Berner Glove	7.5	sterile	50pcs/pack
GL3025	Berner Glove	8.0	sterile	50pcs/pack
GL3030	Berner Glove	8.5	sterile	50pcs/pack
GL3040	Berner Glove	9.0	sterile	50pcs/pack

For more information or to order contact our expert sales team free on
0800 0328 428

or visit our website www.helapet.co.uk