



# PROFILE 3000™

## STERILE NITRILE SURGICAL GLOVES

PNHS



**Profile 3000™** sterile powderfree nitrile surgical gloves are the first nitrile surgical glove to conform to all CE Medical Gloves Standards. The gloves are hand specific, anatomically shaped for comfort and textured for excellent grip. They have superior sensitivity and the beaded cuff ensures that the gloves have a secure fit at the wrist.

Our nitrile surgical gloves are 100% latex free. This totally removes the possibility of allergic reactions from natural rubber proteins for both users and patients. They are substantially stronger than latex surgical gloves and are resistant to a wide range of chemicals.

These gloves are 300mm (12") in length for maximum protection. They have excellent tactility and higher puncture resistance than latex gloves. They are available in a wide range of sizes to suit all users.

*"Comfort, Strength & Dexterity!"*

**HELAPET**

Helapet Limited | Circle Business Centre | Blackburn Road | Houghton Regis | Bedfordshire | LU5 5DD | United Kingdom  
Telephone +44 (0) 1582 501980 | Fax +44 (0) 1582 501981 | e-mail sales@helapet.co.uk | www.helapet.co.uk

Helapet is a CliniMed® Group Company



# PROFILE 3000™

## STERILE NITRILE SURGICAL GLOVES

### PNHS

<b>Colour:</b>	White
<b>Length:</b>	300mm (12")
<b>Material:</b>	Nitrile (Synthetic Latex)
<b>Type:</b>	Hand specific
<b>Surface:</b>	Textured
<b>Other Features:</b>	Powder-free
<b>Protein Level:</b>	Nitrile contains no natural rubber proteins
<b>Sterilisation Method:</b>	Gamma Irradiation
<b>Packaging:</b>	Boxed;
<b>AQL:</b>	1.5 for holes (1 000ml water test)
<b>Physical Properties:</b>	Manufactured to exceed the requirement on the European Medical Gloves Standards EN455 (Parts 1, 2 & 3)
<b>Quality Systems:</b>	Manufactured in a facility holding the following registrations: ISO 9001:2000 and Directive 93/42/EEC Annex V
<b>Sampling &amp; Inspection:</b>	In accordance with ISO 2859
<b>Storage Instructions:</b>	Store in a cool dry place away from direct sunlight and heat (5 to 35)
<b>Shelf Life:</b>	Five years from date of manufacture
<b>Bar Codes:</b>	Each size carries a unique bar code on the inner box (EAN 13) and the carton (ITF 14) to allow automatic entry into Stock Control Systems
<b>Lot Number:</b>	This appears on both the sealed pouch, the inner dispenser box and the outer carton to provide traceability
<b>Packaging:</b>	One cuffed pair per inner wrapper; One inner wrapper per sealed pouch; 50 pairs per dispenser box; 4 dispenser boxes per outer carton (200 pairs)
<b>Regulatory Classification:</b>	Class IIA
<b>Product Codes:</b>	Each size is given an individual code for order purposes. It follows:

Size:	Order code:
6.5	GN9465
7	GN9470
7.5	GN9475
8	GN9480
8.5	GN9485
9	GN9490

*"Comfort, Strength & Dexterity!"*