



# BioClean Synergy™

## NON-STERILE AMBIDEXTROUS NITRILE GLOVES



**BioClean Synergy™** non-sterile powder-free ambidextrous nitrile gloves are one of the most technically advanced gloves available today. They are packed 100 pieces per double bag, to meet the needs of cleanroom users, and are suitable for use in non-sterile ISO Class 4 environments.

Designed with a smooth surface which makes them ideal as an under-glove where double-donning is an important factor, they are suitable for use on their own where extra sensitivity is required.

They are latex free, thus eliminating Type 1 latex allergies. They are also accelerator and sulphur free, thereby virtually eliminating Type 4 allergies.

These ambidextrous nitrile gloves have beaded cuff for added strength, easier donning and secure fit at the arm and are 300mm (12") in length for added protection in cuff area.

Available in a range of sizes to suit all users, they also comply with the European PPE Category 3 Complex Design.





# BioClean Synergy™

## NON-STERILE AMBIDEXTROUS NITRILE GLOVES

Colour:	Green.
Length:	300 mm (12").
Material:	Nitrile.
Type:	Ambidextrous.
Surface:	Microtextured.
Other Features:	Powder-free, latex free beaded cuff, double donnable.
Protein Level:	Zero, Nitrile contains no natural latex proteins.
Packing:	Bagged.
AQL:	1.5 for holes (1000ml water test).
Particulates:	Processed in a cleanroom environment to ensure compatibility with ISO Class 4 operating conditions.
Extractables:	Extended leaching in pure de-ionised water ensures low levels of chlorides and extractable metals.
Physical Properties	Manufactured to exceed the requirements of the PPE Gloves Standard EN420:2003 and the European Medical Gloves Standard EN455 (Parts 1,2 & 3).
Quality Systems:	Manufactured in a facility holding ISO 9001:2000 and 89/686/EEC Article 11b certifications.
Sampling & Inspection:	In accordance with ISO 2859.
Storage Instructions:	Store in a cool, dry place (5°C - 35°C) away from direct sunlight and florescent light.
Shelf Life:	Five (5) years from date of manufacture.
Bar Codes:	Each size carries a unique bar code on the outer bag (EAN 13) and carton (ITF 14).
Lot Number:	This appears on the sealed outer bag and carton to provide traceability.
Packaging:	100 pieces per sealed inner PE bag; One sealed inner PE bag per sealed outer PE bag; 10 outer PE bags per carton liner; One Carton liner per carton (1000 pieces).
Compatibility:	ISO Class 4 Compatible.
Regulatory Classification:	Category 3 - Complex Design Personal Protective Equipment conforming to: EN 420:2003 EN 374-1,2 & 3:2003 EN 388:2003



\	8	7	0
GU2060	XS	Hand size 6	290mm
GU2065	S	Hand size 6.5	290mm
GU2075	M	Hand size 7.5	290mm
GU2085	L	Hand size 8.5	290mm
GU2090	XL	Hand size 9	290mm
GU2095	XXL	Hand size 9.5	290mm

*"The Underglove with Great Tactility!"*



GB04/61092