



Product Specification

Product Name:	BioClean Emerald
Product Description:	Sterile Cleanroom Gloves
Product Code:	BENS 6.0 – 9.0 (in half sizes) plus 10.0

Glove Design:

Material	Nitrile	Colour	Green
Former	Hand-specific	Powder	Powder free

Spec \ Size	6.0	6.5	7.0	7.5	8.0	8.5	9.0	10.0
Total length, mm (min)	290	290	290	290	290	290	290	290

Palm width, mm (± 5mm)	79	83	89	95	103	110	112	120
------------------------	----	----	----	----	-----	-----	-----	-----

Single wall Thickness, mm (min)	Cuff	Palm	Finger Tip
	0.04	0.06	0.08

Protein Content (µg/g)	Zero	Powder Content (mg/glove)	2 (max)
------------------------	------	---------------------------	---------

Physical Properties		
Reference Standard	Before Ageing	After Ageing
	Force at break (N)	Force at break (N)
EN 455	>3.6	>3.6

Quality Assurance Standard						
Sampling Plan Used Single Normal						
Random Sampling based on ISO 2859	Characteristic	Dimension	Physical Properties	Water Tight	Visual Defect (Major)	Visual Defect (Minor)
	Inspection Level	S-2	S-2	G-1	G-1	G-1
	AQL	4.0	4.0	1.5	2.5	4.0

Packaging Instructions:

- 1 pair per inner wallet per sealed easy tear pouch
- 10 pairs per sealed inner bag
- 20 inner bags per carton liner
- 1 carton liner per carton (200 pairs)

Contact Details:

Europe	North America	Asia
NITRITEX LTD Minton Enterprise Park, Oaks Drive, Newmarket, Suffolk, CB8 7YY United Kingdom Tel: +44 (0) 1638 663338 Fax: +44 (0) 1638 668890 Email: sales@nitritex.com	NITRITEX CANADA LTD 7030 Woodbine Avenue, Suite 500, Markham, Ontario L3R 6G2 Canada Tel: +1 (905) 946 9539 Fax: +1 (905) 946 8584 Email: nitritex@cs.com	NITRITEX ASIA PTE LTD 11 Keng Cheow Street, #03-09 Riverside Piazza, Singapore 059608, Republic of Singapore Tel: +65 6536 9123 Fax: +65 6438 5159 Email: nitritex@pacific.net.sg