

SteriClean®

BIO+ II
biocidal disinfectant



together we're stronger

- Double bagged and sterile
- Active blend of mix of chlorhexidine and quaternary ammonium compounds
- Sporistatic to vegetative fungal spores
- Compatible with **BIO+ I** and SteriClean® alcohols
- Broad spectrum activity

the NEW SteriClean®
cleaning regime from
HELAPET

Clean smarter with the NEW SteriClean® rotational regime from Helapet

Introducing BIO+ II

SteriClean® introduce BIO+ II, to compliment the existing SteriClean® denatured ethanol and isopropyl alcohol sprays. Used in conjunction with SteriClean® alcohol, wipes and neutral detergents, this helps create the ideal combination for maximum effect against microbial contamination commonly found in cleanroom environments.

BIO+ II contains an active mix of chlorhexidine, and quaternary ammonium compounds, blended with WFI filtered to 0.2 micron and packed in a class 10,000 cleanroom. Supplied double bagged and sterile for aseptic areas, SteriClean® BIO+ II is terminally sterilised by gamma irradiation at no less than 25kgy.

- Ready to use straight from the dispenser and validated in use shelf life of 3 months.
- Non-irritant and non-corrosive with jet or spray for wider coverage.

Test results

BS EN 1276 (1997)

E-Coli, E Hirae, P. Aeruginosa, S aureus
Sporistatic to fungal spores

For more information and to try SteriClean® cleaning in your unit, please call our sales team on

0800 0328 428

or visit www.helapet.co.uk



Order number

XX6602

Size

sterile 0.9L 6 bottles/case

HELAPET