

# HELAPET

Filter Devices • Medical & Cleanroom Disposables

**FOR USE WITH  
ISOLATORS & CABINETS**



## **Berner Chemotherapy Gloves**

ANATOMICAL • STERILE • POWDER FREE

GOOD SENSITIVITY • NO ALLERGENS DETECTED

**Freefone Orders: 0800 0 328 428**

### **Helapet Limited**

Circle Business Centre  
Blackburn Road  
Houghton Regis  
Bedfordshire  
LU5 5DD  
England

**Tel: 01582 501980**

**Fax: 01582 501981**

**www.helapet.co.uk**

**E-mail: info@helapet.co.uk**

A CliniMed® Group Company



CERTIFICATE No. FS 52480  
BS EN ISO 9002  
BS EN ISO 13485

# Disposable products for working with CytoToxics

## Berner Gloves for Chemotherapy Preparation

This protective glove has been developed for handling chemotherapeutic agents currently used in the course of medical treatment. In view of the wide range of chemotherapeutic agents on the market Berner cannot give a universal guarantee of resistance or protection. When using Carmustin or Thiotepa for example the gloves should be changed after every usage. A visible contamination demands an immediate change of the gloves. These gloves contain materials, which are known to possibly cause allergic reactions: Latex, Carbamate. Physical requirement tests are specified in EN 374 and EN 388: abrasion resistance at least 100 cycles, cut resistance index 1.04.

Permeation By Cytotoxics : Carmustin class 4 (> 120 min), 5-fluorouracil (> 360 min).

Product	Code	Qty
Latex, Powder Free, anatomical, size 6.5	GL3010	50 prs/pack
Latex, Powder Free, anatomical, size 7	GL3015	50 prs/pack
Latex, Powder Free, anatomical, size 7.5	GL3020	50 prs/pack
Latex, Powder Free, anatomical, size 8	GL3025	50 prs/pack
Latex, Powder Free, anatomical, size 8.5	GL3030	50 prs/pack
Latex, Powder Free, anatomical, size 9	GL3040	50 prs/pack

## Nitrile Gloves for latex sensitive Chemotherapy Administration

Non latex, powder free sterile anatomical gloves that have high resistance to all commonly used chemotherapeutic agents (data available on request).

Product	Code	Qty
Nitrile, anatomical, powder free, 12", size 6.5	GN9465	50 prs/pack
Nitrile, anatomical, powder free, 12", size 7	GN9470	50 prs/pack
Nitrile, anatomical, powder free, 12", size 7.5	GN9475	50 prs/pack
Nitrile, anatomical, powder free, 12", size 8	GN9480	50 prs/pack
Nitrile, anatomical, powder free, 12", size 8.5	GN9485	50 prs/pack
Nitrile, anatomical, powder free, 12", size 9	GN9490	50 prs/pack

## Gowns

Berner gown, polymer coated breathable, non-woven Polypropylene fabric with knitted cuff and comfortable collarless neck with easy to use velcro fastening. Small, medium & large sizes

Product	Code	Qty
Gown, Polymer Coated, Elasticated Cuff, Small	AB2010	10 pcs/pack
Gown, Polymer Coated, Elasticated Cuff, Medium	AB2020	10 pcs/pack
Gown, Polymer Coated, Elasticated Cuff, Large	AB2030	10 pcs/pack

## Oversleeves

Berner protective oversleeves with low lint polymer coating & knitted cuff, sterile. Cone-style cut ensures wearer comfort. Offering optimum personal and product protection. Liquid impermeable coating on the arms. Made from low linting materials for clean environment safety.

Product	Code	Qty
Oversleeve, Polymer coated, standard size	AB0001	40 prs/pack

## Workmats

Berner workmats: To contain spilt liquids and protect the work environment.

Lint free 3 ply absorbent mat (blue). Liquid permeable top layer, highly absorbent intermediate layer and a protective liquid impermeable 1 mil antistatic Polyethylene base layer.

Helapet workmats: Cellulose 1 ply mat with polythene backing, white. Standard and custom cut sizes available.

Product	Code	Qty
Cytotoxic Workmat, 3 ply, 26 x 42cm	WP8001	50 pcs/pack
Cytotoxic Workmat, 3 ply, 41 x 56cm	WP8101	40 pcs/pack
Sterile Pharmacy absorbent workmat, 20 x 25cm	XX2025	200 pcs/pack
Sterile Pharmacy absorbent workmat, 28 x 46cm	XX2846	100 pcs/pack
Pharmacy absorbent workmat, 46 x 57cm	WP4657	50 pcs/pack