

Visimark

A Complete Surgical Skin Marking Kit With A *Measurable* Difference



Peel: the sterile “see through” fibre free pouch.

Measure: with the convenient, flexible ruler included (metric & imperial)

Mark: with VISIMARK, allowing you to draw from a fine line to a medium line in scrub resistant ink.

Also available in an

- ◆ ultra fine tip
- ◆ broad tip



CE marked, sterile

For best results: Remove all oil from skin before marking, using alcohol. Work a heavy line mark into the skin using several strokes of the pen.

Product	Code	Qty
Vismark Surgical Skin Marker Ultrafine	PT1402	50 pcs/pack
Vismark Surgical Skin Marker Regular	PT1400	50 pcs/pack
Vismark Surgical Skin Marker Broad	PT1404	50 pcs/pack

Freefone Orders: 0800 0 328 428

Helapet Limited
 Circle Business Centre
 Blackburn Road
 Houghton Regis
 Bedfordshire
 LU5 5DD
 England

Tel: 01582 501980
Fax: 01582 501981
www.helapet.co.uk
E-mail: info@helapet.co.uk

A CliniMed® Group Company



CERTIFICATE No. FS 52480
 BS EN ISO 9002
 BS EN ISO 13485