

INSTANT NO CONTACT TEMPERATURE READING



Helapet is proud to announce the launch in the UK of a professional clinical thermometer that takes temperature from a distance instantly, and without contact with the patient. Designed and manufactured by the French Company Visiomed, the ThermoFlash® PRO LX261 has been developed using the latest infrared technology. This technology allows temporal artery temperature to be taken at a distance of about five centimetres away from the forehead. As there is no contact with the patient, there is no risk of cross contamination and is ideal for areas where infection prevention is paramount.

ThermoFlash® is easy to use and extremely accurate. It has been clinically tested in France and Helapet has had an independent validation study commissioned to verify this. The revolutionary concept allows you to capture, analyse and display a person's temperature based on the energy that person releases. All this is obtained in a fraction of a second.

This product will save you time and money. It uses no consumable probes, making it ideal for taking the temperature of newborn babies, young children, the elderly and for critical care patients. Because it can be used, even when asleep, it minimises disturbance and discomfort for patients.

It is a fully CE marked, Class 2a medical device and comes with a two year guarantee. It has a fully rechargeable battery which can carry out over ten thousand tests per charge. The battery life is over 3.6 million cycles.

ThermoFlash® is a versatile device. Not only does it record body temperature, but it can also monitor room temperature and that of liquids, such as bath water and babies' feed bottles.

For more information or to request a demonstration, please call 0800 0328 428 or visit our website at www.helapet.co.uk

Helapet will be demonstrating ThermoFlash® at **Infection Prevention 2011** at the Bournemouth International Centre, 19th-21st September. Be sure to come and visit us on Stand **W3**.