

STOP PRESS ** STOP PRESS **** STOP PRESS **** STOP PRESS**

Company: Helapet Ltd
Contact: Ruth Hewitt
Telephone: 01582 501983
Date: 16th February 2010

Title: TimestripPlus™ for TPN

New TimestripPlus™ is a cost effective method of making sure that Total Parenteral Nutrition (TPN) made up in pharmacy is delivered to the patient with no break in the cold chain.

There is currently no other inexpensive method of tracking the temperature of cold chain items once they have left the pharmacy or validated refrigeration unit.

Total Parenteral Nutrition (TPN) is a typical example of a product manufactured in aseptic pharmacy that must be kept in a validated refrigerated environment and not left at room temperature for lengthy periods of time which could result in compromising the safety of immune deficient patients.

Simply activate the TimestripPlus™ at room temperature and apply to the TPN bag using the self adhesive backing. The TPN bag should be kept below 8 deg C (refrigerated) until it is required.

When the TPN bag is delivered to the patient and reaches 12 deg C, the indicator will travel along the window displaying the amount of time that the bag of TPN has reached 12 deg C. If the bag is unused and returned to the fridge (8 deg C) the indicator will stop.

Using TimestripPlus™ for TPN offers peace of mind that the bag has not had a temperature breach before use and significant cost savings can be achieved if unused and chilled product can be returned to pharmacy for re-issue.

TimestripPlus™ is the only device that not only records a breach in the cold chain but indicates how long that breach occurred.

For more information and FREE samples please call **0800 0328 428** or visit **www.helapet.co.uk**

Ends
273 words